


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4'-Iodoacetophenone - $\geq 97\%$

Product Number: **357804**
CAS Number: 13329-40-3
MDL: MFCD00045320
Formula: C₈H₇IO
Formula Weight: 246.05 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White to Yellow and White to Orange and White to Brown
Appearance (Form)	Powder or Crystal or Chunk(s)
Infrared spectrum	Conforms to Structure
Titration by AgNO ₃ After Oxygen Combustion	97.0 - 103.0 %
Purity (GC)	$\geq 97.0 \%$

Specification: PRD.0.ZQ5.10000079618

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.