


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Sodium 5-methylpyridine-2-sulfinate - ≥95%

Product Number: 900790



Formula: C<sub>6</sub>H<sub>6</sub>NNaO<sub>2</sub>S  
Formula Weight: 179.17 g/mol

TEST	Specification
Appearance (Color)	White to Beige
Appearance (Form)	Powder or Crystals
<sup>1</sup> H NMR Spectrum	Conforms to Structure
Carbon Content	38.2 - 42.3 %
Nitrogen Content	7.4 - 8.3 %
Sulfur Content	16.9 - 18.8 %
Water (by Karl Fischer)	
Specification: Report Result	
Purity (HPLC)	≥ 95.0 %
ICP Major Analysis	Confirmed
Confirms Sodium Component	

Specification: PRD.1.ZQ5.10000084253

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.