


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Product Number: 159220
CAS Number: 1609-47-8
MDL: MFCD00009106
Formula: C₆H₁₀O₅
Formula Weight: 162.14 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Proton NMR spectra	Conforms to Structure
Purity (Titration by HCL04)	96.5 - 103.5 %
Color Test	≤ 20 APHA

Specification: PRD.0.ZQ5.10000028191

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.