

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Trifluoromethanesulfonic acid - reagent grade, 98%

Product Number:

158534

CAS Number:

1493-13-6

MDL:


MFCD00007514

Formula:

CHF₃O₃S

Formula Weight:

150.08 g/mol

**TEST****Specification**

Appearance (Color)

Colorless to Yellow

Appearance (Form)

Liquid

Titration with NaOH

97.5 - 102.5 %

Color Test

≤ 400 APHA

Specification: PRD.2.ZQ5.10000004552

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.