

3050 Spruce Street, Saint Louis, MO 63103, USA

ZrO 2 / Y2O3

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Zirconium(IV) oxide-yttria stabilized - submicron powder, 99.9% trace metals basis (purity excludes ~2% HfO2)

Product Number: 464228

CAS Number: 114168-16-0 MDL: MFCD00145546

MDL: MFCD001 Formula: O5Y2Zr

Formula Weight: 349.03 g/mol

TEST Specification

Appearance (Color) White

Appearance (Form) Powder

X-Ray Diffraction Conforms to Structure

d50 (Microtrac X100)

Size

Specification: Report Result

Purity Meets Requirements

99.9% Based on Trace Metal Analysis Excluding Hafnium

Trace Metal Analysis ≤ 1500.0 ppm

Specification: PRD.2.ZQ5.10000031753

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.