

3050 Spruce Street, Saint Louis, MO 63103, USA

GeCl<sub>4</sub>

Meets Requirements

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

## **Product Specification**

Product Name:

Germanium(IV) chloride - 99.99% trace metals basis

 Product Number:
 208450

 CAS Number:
 10038-98-9

 MDL:
 MFCD00011029

Formula: Cl4Ge

Formula Weight: 214.45 g/mol

TEST Specification

Appearance (Color) Colorless

Appearance (Form) Liquid

Purity 99.99% Based On Trace Metals Analysis

Trace Metal Analysis < 150.0 ppm

Specification: PRD.1.ZQ5.10000025609

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.