

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)

Email USA: [techserv@sial.com](mailto:techserv@sial.com)

Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)


## Product Specification

Product Name:  
Naphthol AS-BI Phosphoric acid solution

Product Number: **3871**

MDL: MFCD00004058

Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Colorless to Pale Yellow
Appearance (Form)	Solution
Appearance (Turbidity)	Clear
Essentially free of particles	
Performance	Conforms
Stains within performance characteristics when using Sigma-Aldrich No. 387	
Expiration Date Period	-----
18 months	

Specification: PRD.0.ZQ5.10000065109

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.