

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Cholesteryl linoleate -  $\geq 98\%$  (HPLC; detection at 205 nm)

Product Number:

C0289

CAS Number:

604-33-1

MDL:

MFCD00003644

Formula:


C45H76O2

Formula Weight:

649.08 g/mol

Storage Temperature:

-20 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
100 mg/mL, CHCl <sub>3</sub>	
Infrared spectrum	Conforms to Structure
Purity (HPLC)	$\geq 98\%$

Specification: PRD.0.ZQ5.10000011729

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.