SIGMA-ALDRICH®

sigma-aldrich.com

3050 Spruce Street, Saint Louis, MO 63103, USA Website: www.sigmaaldrich.com Email USA: techserv@sial.com Outside USA: eurtechserv@sial.com

Product Specification

Product Name: EX-CELL* Advanced CHO Fed-batch Medium - liquid

Product Number: 14366C

Storage Temperature: 2 - 8 °C

TEST	Specification	
Appearance (Turbidity)	Clear	
Appearance (Form)	Solution	
with no precipitates		
Appearance (Color)	Report Only (Color)	
рН	7.2 - 7.4	
Osmolality	280 - 320 mOs/kg	
Sterility by USP guidelines	Sterile	
Endotoxin Level		
Record		
Growth Promotion (cells/ml)		
Record		
Viability		
Record		
Site of Manufacturing	Manufacturing Site	

Remarks:

Date of manufacture(DDMMMYY):

Specification: PRD.4.ZQ5.10000057996

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.