

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Anti-Leukemia Inhibitory Factor antibody produced in goat - IgG fraction of antiserum, lyophilized powder

Product Number: **L9152**

Storage Temperature: -20 °C

TEST	Specification
Formulation Lyophilized from Phosphate Buffered Saline, containing 5% Trehalose.	Conforms
Bioactivity (LIF Neutralization, DA-1A Cells) ND50: 80-480 ng/mL	Conforms
Working Dilution by Direct Immunoblotting 0.1 ug/mL will allow detection of recombinant mouse LIF	Conforms
Endotoxin Assay < 0.1 EU/ug	Conforms

Specification: PRD.1.ZQ5.10000049770

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.