


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Chloro-5-iodo-3-nitropyridine - 97%

Product Number: **720054**
CAS Number: 426463-05-0
Formula: C₅H₂ClIN₂O₂
Formula Weight: 284.44 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White to Yellow
Appearance (Form)	Powder to Crystalline Powder
Infrared spectrum	Conforms to Structure
Titration by AgNO ₃ After Oxygen Combustion	96.5 - 103.5 %
Purity (GC)	≥ 96.5 %

Specification: PRD.0.ZQ5.10000028160

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.