

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

4-Chloro-2-methylbenzyl alcohol - 97%

Product Number:

647551

CAS Number:

129716-11-6

MDL:


MFCD06201140

Formula:

C₈H₉ClO

Formula Weight:

156.61 g/mol



TEST

Specification

Appearance (Color)

White to Off-White

Appearance (Form)

Powder

Infrared spectrum

Conforms to Structure

Titration by AgNO₃ (CLCO)

96.5 - 103.5 %

After Oxygen Combustion

Purity (GC)

≥ 96.5 %

Specification: PRD.0.ZQ5.10000045246

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.