

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Carbon nanotube, multi-walled - as-produced cathode deposit, >7.5% MWCNT basis, O.D. × L 7–15 nm × 0.5–10 μm

Product Number:

412988

CAS Number:

308068-56-6



TEST	Specification
------	---------------

Appearance (Color)	Black
Appearance (Form)	Conforms to Requirements
Cylinders or Powdery Solid	
Graphitic Layers	Confirmed
Approximately 5-50 Graphitic Layers	
Note	Confirmed
Product Identity Assured by Supplier Based On Manufacturing Method	

Specification: PRD.0.ZQ5.10000087145

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.