

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)

Email USA: [techserv@sial.com](mailto:techserv@sial.com)

Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

Product Specification

Product Name:  
Hanahan's Broth (SOB Medium) – powder microbial growth medium

Product Number: H8032

TEST	Specification
Appearance (Color)	Off-White to Yellow
Appearance (Form)	Powder
Solubility	Conforms
Clear, Amber Solution; 2.8% w/v H2O	
pH	6.8 - 7.2
2.8% w/v H2O	
Cultural Response	Pass
as described in data sheet	

Specification: PRD.1.ZQ5.10000021188

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.