


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Cholesteryl erucate - $\geq 95\%$

Product Number: C7259
CAS Number: 24516-39-0
MDL: MFCD00037704
Formula: C₄₉H₈₆O₂
Formula Weight: 707.21 g/mol
Storage Temperature: -20 °C



TEST	Specification
------	---------------

Appearance (Color)	White
Appearance (Form)	Powder
Purity (GC)	$\geq 95\%$
Specific Rotation	-22 - -19 °
C = 1.0 in CHCl ₃ at 25 deg C.	
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
50 mg/mL, CHCl ₃	

Specification: PRD.0.ZQ5.10000094618

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.