

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
PALLADIUM ON CARBON, EVONIK NOBLYST(R) P&

Product Number: **908061**

TEST	Specification
Appearance (Color)	Black
Appearance (Form)	Powder
ICP Major Analysis Confirms Palladium Component	Confirmed
Loss on Drying Calculated to dry basis	45.0 - 65.0 % wt
Palladium (Pd) Calculated to dry basis	4.5 - 5.5 % wt
Activity Cinnamic acid hydrogenation ml/min.g	≥ 575
Product of Supplier Evonik	Confirmed
Registered Trademark Noblyst is a registered trademark of Evonik Degussa GmbH	Note
Recommended Retest Period 2 Years	-----

Specification: PRD.0.ZQ5.10000100341

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.