

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

1,3-Dihydro-1,3,3-trimethylspiro[2H-indole-2,3'-[3H]naphth[2,1-b][1,4]oxazine] - ≥98.0% (HPLC)

Product Number:

322547

CAS Number:

27333-47-7

MDL:


MFCD00013318

Formula:

C₂₂H₂₀N₂O

Formula Weight:

328.41 g/mol

**TEST****Specification**

Appearance (Color)

Yellow to Tan

Appearance (Form)

Conforms to Requirements

Powder or Crystalline Powder

Infrared spectrum

Conforms to Structure

Purity (HPLC)

≥ 98 %

Specification: PRD.0.ZQ5.10000054837

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.