

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Lipopolysaccharides from Escherichia coli 0111:B4 - FITC conjugate

Product Number:

F3665

TEST	Specification
Appearance (Form)	Lyophilized Powder
Appearance (Colour)	Yellow to Orange
Solubility (Solvent)	Water
Solubility (Conc)	4.90 - 5.10 mg/ml
Solubility (Turbidity)	Faint Hazy to Hazy
Solubility (Color)	Pale Yellow to Yellow
Potency (Sample EU/mg)	≥ 100000 EU/mg
Conj. (Ratio FITC/LPS ug/mg)	2.00 - 10.00 µG/mg
Recommended Retest Period	-----
2 Years	

Specification: PRE.0.ZQ5.10000001383

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.