


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
R-MOP

Product Number: 693537
CAS Number: 145964-33-6
MDL: MFCD00269686
Formula: C₃₃H₂₅OP
Formula Weight: 468.52 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
Solid, Powder, Crystals, Crystalline Powder and/or Chunks
Infrared spectrum
Phosphorus NMR Spectra
Optical Rotation
c = 0.5%; Chloroform
Purity (HPLC)
Note
Product of Takasago International Corporation

White to Beige
Conforms to Requirements

Conforms to Structure
Conforms to Structure
89.0 - 99.0 deg

≥ 94.0 %
Note

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.