

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

VIRAL RNA EXTRACTION BUFFER

Product Number:

VRE100

Storage Temperature:

-20 °C

TEST	Specification
Appearance Turbidity	Clear
Appearance (Color)	Colorless to Light Yellow
Appearance (Form)	Liquid
Suitability	Pass
Suitable for use in qRT-PCR reactions	
Suitability-SARS-CoV-2	None Detected
No amplification in 45 cycles	

Specification: PRD.1.ZQ5.10000103457

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

