

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

D-(–)-3-Phosphoglyceric acid disodium salt – ≥93%, powder

Product Number: P8877
CAS Number: 80731-10-8
MDL: MFCD00069569
Formula: C₃H₅Na₂O₇P
Formula Weight: 230.02 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear
50 mg/mL, H ₂ O	
Water (by Karl Fischer)	≤ 14 %
Sodium (dry basis)	18 - 22 %
By ICP	
Enzymatic Purity (Dry Basis)	≥ 93 %
Solvent Content (by NMR)	≤ 6 %

Specification: PRD.1.ZQ5.10000018424

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.