

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification


Product Name:
Acetone - Laboratory Reagent, ≥99.5%

Product Number: 179973

CAS Number: 67-64-1

Formula: C₃H₆O

Formula Weight: 58.08 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Refractive index at 20 ° C	1.357 - 1.361
Purity (GC)	≥ 99.50 %
Water (by Karl Fischer)	≤ 0.5 %
Residue on Evaporation	≤ 0.002 %

Specification: PRD.1.ZQ5.10000018591

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.