

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium cyclohexyltrifluoroborate - 95%

Product Number:

723576

CAS Number:


446065-11-8

Formula:

C₆H₁₁BF₃K

Formula Weight:

190.06 g/mol



TEST

Specification

Appearance (Color)

White to Faint Yellow and Faint Brown

Appearance (Form)

Powder or Crystals

Proton NMR spectrum

Conforms to Structure

Fluorine NMR Spectra

Conforms to Structure

Carbon

35.6 - 40.2 %

ICP Major Analysis

Confirmed

Confirms Boron and Potassium Components

Specification: PRD.1.ZQ5.1000022488

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.