

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Palladium(II) chloride - ≥99.9%

Product Number: 520659
CAS Number: 7647-10-1
MDL: MFCD00003558
Formula: Cl₂Pd
Formula Weight: 177.33 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
Orange-Brown or Red-Brown	
Appearance (Form)	Powder
X-Ray Diffraction	Conforms to Structure
Palladium (Pd)	59.0 - 60.0 %
Purity	Meets Requirements
> = 99.9% Based on Trace Metals Analysis (Excluding Sodium)	
Trace Metal Analysis	≤ 1000.0 ppm

Specification: PRD.1.ZQ5.10000049685

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.