

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium metatitanate – –200 mesh

Product Number: 401307
CAS Number: 12034-36-5
MDL: MFCD00014243
Formula: Na₂O₇Ti₃
Formula Weight: 301.58 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
ICP Major Analysis	Confirmed
Confirms Na and Ti Components	
Size	≤ 75.0 micron
Average Particle Size	

Specification: PRD.0.ZQ5.10000021934

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.