

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Anti-Heme Oxygenase-1 antibody produced in rabbit - IgG fraction of antiserum, buffered aqueous solution

Product Number: **H4535**

MDL: MFCD01092544

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Formulation	Conforms
Solution in phosphate buffered saline, containing 30% glycerol, 0.5% BSA, and 0.01% thimerosal.	
Concentration	≥ 0.5 mg/ml

Specification: PRD.0.ZQ5.10000056967

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.