

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

3-Aminopyrrolidine dihydrochloride - 98%

Product Number:

404624

CAS Number:


103831-11-4

MDL:

MFCD00060176

Formula Weight:

159.06 g/mol



### TEST

### Specification

Appearance (Color)

White to Light Yellow and Light Beige to Beige

Appearance (Form)

Powder or Crystals or Chunk(s)

Infrared Spectrum

Conforms to Structure

Carbon Content

29.4 - 31.0 %

Nitrogen Content

17.1 - 18.1 %

Specification: PRD.1.ZQ5.10000103195

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

