


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1,3,5-Triallyl-1,3,5-triazine-2,4,6(1H,3H,5H)-trione - 98%

Product Number: 114235
CAS Number: 1025-15-6
MDL: MFCD00006554
Formula: C₁₂H₁₅N₃O₃
Formula Weight: 249.27 g/mol



TEST

Specification

Appearance (Color)	Colorless or White
Appearance (Form)	Conforms to Requirements
Viscous Liquid or Solid	
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 97.5 %

Specification: PRD.0.ZQ5.10000018413

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.