

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Methyl methacrylate - 99%, contains ≤ 30 ppm monomethyl ether hydroquinone as inhibitor

Product Number:

M55909

CAS Number:

80-62-6

MDL:

MFCD00008587

Formula:


C₅H₈O₂

Formula Weight:

100.12 g/mol

Storage Temperature:

2 - 8 °C

**TEST****Specification**

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Refractive index at 20 °C

1.412 - 1.416

Infrared spectrum

Conforms to Structure

Purity (GC)

 ≥ 98.5 %

MEHQ as Inhibitor

8 - 30 ppm

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.