

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Protease Inhibitor Cocktail – for use in tissue culture media, DMSO solution

Product Number: P1860

Storage Temperature: -20 °C

TEST	Specification
Appearance (Form)	Clear Solution
Appearance (Colour)	Colorless to Faint Yellow
Suitability (Result)	Pass
Vial Content (Size)	1 - 1
Vial Content (Unit Of Measure)	ml
Vial Content (Actual Content)	≥ 1.00
Recommended Retest Period	-----
4 Years	

Specification: PRE.0.ZQ5.10000001989

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.