

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
TEMPO - purified by sublimation, 99%

Product Number: 426369
CAS Number: 2564-83-2
MDL: MFCD00009599
Formula: C₉H₁₈NO
Formula Weight: 156.25 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)	Orange to Red
Appearance (Form)	Solid or Powder or Crystals or Chunks
Infrared spectrum	Conforms to Structure
Melting Point	
Purity (GC)	≥ 98.5 %
Solubility (Turbidity)	Clear
50mg/ml (5%); Ethanol	
Solubility (Color)	Orange to Red
Recommended Retest Period	-----
2 Years	

Specification: PRD.0.ZQ5.10000048395

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.