

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-(4-Ethyl-2-pyridyl)-1H-imidazole - 98%

Product Number: 757195

Formula: C₁₀H₁₁N₃
Formula Weight: 173.21 g/mol

TEST	Specification
------	---------------

Appearance (Color)	White to Faint Beige
Appearance (Form)	Powder or Solid or Crystals or Crystalline Chunk(s) or Chunk(s)
Infrared spectrum	Conforms to Structure
Carbon	67.6 - 71.1 %
Nitrogen	23.6 - 24.9 %
Purity (GC)	≥ 97.5 %

Specification: PRD.0.ZQ5.10000028450

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.