

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name: CAPS - ≥99%

 Product Number:
 C2632

 CAS Number:
 1135-40-6

 MDL:
 MFCD00003837

 Formula:
 C9H19NO3S

 Formula Weight:
 221.32 g/mol

TEST

Appearance (Color)
Appearance (Form)
Solubility (Color)
Solubility (Turbidity)

at 15 g + 85 mL, H2O

Infrared spectrum
Water (by Karl Fischer)
Carbon (anhydrous)
Nitrogen (anhydrous)
A290 UV absorbance

15% W/W ICP-MS

(Total: Ag, As, Bi, Cd, Cu, Hg, Mo, Pb, Sb, Sn)

Titration with NaOH

(Anhydrous)

Recommended Retest Period

3 years

Specification: PRD.3.ZQ5.10000052001

Specification

White

Crystalline Powder

Colorless Clear

Conforms to Structure

< 1 %

48.3 - 49.3 %

6.1 - 6.5 %

< 0.05

< 5 ppm

> 99 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.