


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Resveratrol - ≥99% (HPLC)

Product Number: R5010
CAS Number: 501-36-0
MDL: MFCD00133799
Formula: C₁₄H₁₂O₃
Formula Weight: 228.24 g/mol
Storage Temperature: -20 °C



TEST

Specification

Appearance (Color)	White to Off-White
Appearance (FORM)	Powder
Solubility (Color)	Faint Yellow
Solubility (Turbidity)	Clear
50 mg/ml, Acetone	
Infrared spectrum	Conforms to Structure
Carbon	73.0 - 74.3 %
Purity (HPLC)	≥ 99 %

Specification: PRD.2.ZQ5.10000010053

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.