


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Magnolol - $\geq 95\%$ (HPLC), from plant

Product Number: **M3445**
CAS Number: 528-43-8



Formula: C₁₈H₁₈O₂
Formula Weight: 266.33 AMU

TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
1.5 mg/ml, MeOH	
Proton NMR spectrum	Conforms to Structure
¹³ C NMR Spectrum	Conforms to Structure
Purity (HPLC)	$\geq 95\%$

Specification: PRE.0.ZQ5.10000020765

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.