

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Sodium DL-lactate solution – syrup, BioReagent, 60 % (w/w), synthetic, suitable for cell culture

Product Number:

L4263

CAS Number:

72-17-3

MDL:

MFCD00065400

Formula:


C₃H₅NaO₃

Formula Weight:

112.06 g/mol

Storage Temperature:

2 - 8 °C



TEST

Specification

Appearance (Color)

Colorless

Appearance (Form)

Viscous Liquid

Appearance (Turbidity)

Clear

Water (by Karl Fischer)

35 - 45 %

Cell Culture Test

Pass

Sodium (Na)

10.2 - 13.7 %

Purity (Enzymatic)

27 - 33 %

L-isomer

Recommended Retest Period

5 years

Specification: PRD.1.ZQ5.10000040909

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.