


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4,5-Dibromo-2-furoic acid - 97%

Product Number: 656291
CAS Number: 2434-03-9
MDL: MFCD00092311
Formula: C₅H₂Br₂O₃
Formula Weight: 269.88 g/mol



TEST	Specification
Appearance (Color)	White to Brown
Appearance (Form)	Crystal to Powder
Infrared spectrum	Conforms to Structure
GC (area %)	≥ 96.5 %

Specification: PRD.1.ZQ5.10000041139

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.