

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Monoclonal Anti-Cytokeratin, pan-FITC antibody produced in mouse - clone C-11, purified immunoglobulin, buffered aqueous solution

Product Number:

F3418

TEST	Specification
------	---------------

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Yellow
IgG Concentration (OD 280nm)	1.5 - 3.0 mg/ml
F/P Molar Ratio	3.0 - 8.0
ICC/IF Cells	PTK2-Marsupial Kidney Cells
ICC/IF Fixation Reagent	Methanol / Acetone (1:1)
ICC/IF-Titer	1:50-1:100
IHC - Tissue	Human Placenta
IHC - Fixation	FFPE-Formalin Fixed Paraffin Embedded
IHC - Titer	1:50-1:100
Recommended Retest Period	-----
3 Years	

Specification: PRE.0.ZQ5.10000016055

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.