

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification


Product Name:
Bis[rhodium($\alpha,\alpha,\alpha',\alpha'$ -tetramethyl-1,3-benzenedipropionic acid)] - 96%

Product Number: **662623**

CAS Number: 819050-89-0

Formula: C₃₂H₄₀O₈Rh₂

Formula Weight: 758.47 g/mol



TEST

Specification

Appearance (Color)

Conforms to Requirements

Green

Appearance (Form)

Conforms to Requirements

Powder, Crystals or Crystalline Powder, Solid and/or Chunks

Infrared spectrum

Conforms to Structure

Carbon

48.4 - 53.0 %

Purity (HPLC)

≥ 95.5 %

Specification: PRD.0.ZQ5.10000036598

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.