


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Coenzyme Q10 - ≥98% (HPLC)

Product Number: **C9538**  
CAS Number: 303-98-0  
MDL: MFCD00042919  
Formula: C<sub>59</sub>H<sub>90</sub>O<sub>4</sub>  
Formula Weight: 863.34 g/mol  
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	Yellow to Dark Orange
Appearance (Form)	Powder
Solubility (Color)	Yellow to Orange
Solubility (Turbidity)	Clear to Slightly Hazy
50 mg/mL, CHCl <sub>3</sub>	
E1%	≥ 162
at 275 nm in EtOH	
Purity (HPLC)	≥ 98 %

Specification: PRD.0.ZQ5.10000048402

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.