


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3,3-Dimethylacrylic acid - 97%

Product Number: D138606
CAS Number: 541-47-9
MDL: MFCD00004366
Formula: C₅H₈O₂
Formula Weight: 100.12 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White to Light Yellow to Light Beige
Appearance (Form)	Powder or Crystals
Infrared Spectrum	Conforms to Structure
Titration with NaOH	96.5 - 103.5 %
GC (area %)	≥ 96.5 %

Specification: PRD.0.ZQ5.10000085104

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.