


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2,2'-Azobis(2-methylpropane) - 97%

Product Number: 452106
CAS Number: 927-83-3
MDL: MFCD00014998
Formula: C₈H₁₈N₂
Formula Weight: 142.24 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)
Appearance (Form)
Proton NMR spectra
Purity
97% Purity By HNMR

Colorless to Yellow
Liquid
Conforms to Structure
Conforms

Specification: PRD.1.ZQ5.10000007634

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.