

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Lipase from Candida sp. – recombinant, expressed in Aspergillus niger

Product Number: L3170

Storage Temperature: 2 - 8 °C

TEST	Specification
Enzyme Activity	≥ 5000
Expressed in LU/g	
Note	-----
Density 1.2 g/ml	

Specification: PRD.0.ZQ5.10000080109

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.