


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
3-Hydroxybutyric acid - 95%

**Product Number:** 166898  
**CAS Number:** 300-85-6  
**MDL:** MFCD00004546  
**Formula:** C<sub>4</sub>H<sub>8</sub>O<sub>3</sub>  
**Formula Weight:** 104.10 g/mol  
**Storage Temperature:** 2 - 8 °C



TEST	Specification
------	---------------

Appearance (Color)	White to Faint Yellow and Colorless to Faint Yellow
Appearance (Form)	Liquid or Solid
Infrared spectrum	Conforms to Structure
Purity (Titration by NaOH)	94.5 - 105.5 %
GC (area %)	≥ 94.5 %

Specification: PRD.2.ZQ5.10000022301

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

