


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Di-tert-butyl N,N-diethylphosphoramidite - technical grade, 93%

Product Number: 365963
CAS Number: 117924-33-1
MDL: MFCD00134307
Formula: C₁₂H₂₈NO₂P
Formula Weight: 249.33 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Appearance (Clarity)	Slightly Hazy to Hazy
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 92.5 %

Specification: PRD.0.ZQ5.10000038826

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.