

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Monoclonal Anti-Na⁺ /K⁺ ATPase (β1 Subunit) antibody produced in mouse - clone M17-P5-F11Product Number: **A278**

MDL: MFCD01633898

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Formulation	Conforms
-------------	----------

Solution in PBS, 1 mg/mL BSA and 0.05% sodium azide

Working Dil. by Western Blot	Conforms
------------------------------	----------

1:1,000 - 1:10,000

Working Dilution	Conforms
------------------	----------

By Immunohistochemistry; 1:100

Specification: PRD.0.ZQ5.10000052212

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.