

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Peroxidase from horseradish – Vetec™ reagent grade

Product Number:

V900503

CAS Number:

9003-99-0

Storage Temperature:

2 - 8 °C

TEST	Specification
------	---------------

Appearance (Colour)

Brown to Very Dark Brown and Red to Very Dark Red and Red-Brown and Very Dark Red-Brown

Appearance (Form)

Powder

Solubility (Colour)

Orange to Dark Orange and Red to Very Dark Red and Orange-Red and Red-Orange and Orange-Brown and Brown-Orange and Dark Red-Orange

Solubility (Turbidity)

Clear

10 mg/mL of 0.1M Phosphate Buffer, pH 6.0

units/mg Solid

≥ 250

One unit will form 1.0 mg Purpurogallin from Pyrogallol in 20 seconds at pH 6.0 at 20 deg C.

Specification: PRE.0.ZQ5.10000011653

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.