

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Titanium nitride - <3 µm

 Product Number:
 595063

 CAS Number:
 25583-20-4

 MDL:
 MFCD00049596

Formula: NTi

Formula Weight: 61.87 g/mol

TiN

TEST

Appearance (Color) Appearance (Form) X-Ray Diffraction Size

Gravimetric Analysis

% Ti

Specification: PRD.0.ZQ5.10000038517

Specification

Brown to Very Dark Brown Powder Conforms to Structure < 3 micron 74.7 - 80.1 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.