


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Iodotrimethylsilane - 97%

Product Number: 195529
CAS Number: 16029-98-4
MDL: MFCD00001028
Formula: C₃H₉Si
Formula Weight: 200.09 g/mol
Storage Temperature: -20 °C

**TEST****Specification**

Appearance (Color)

Conforms to Requirements

Colorless to Yellow to Orange

Appearance (Form)

Liquid

Proton NMR spectra

Conforms to Structure

Purity by Proton NMR

≥ 96.5 %

Specification: PRD.0.ZQ5.10000031161

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.