


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2,4-Diamino-6-hydroxypyrimidine - 96%

Product Number: D19206
CAS Number: 56-06-4
MDL: MFCD00006098
Formula: C₄H₆N₄O
Formula Weight: 126.12 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
Infrared spectrum
Purity (Titration by HCL04)
Solubility (Turbidity)
c = 50 mg/ml (5%); 1N HCl
Solubility (Color)
Purity (HPLC)

White to Beige
Powder or Crystals
Conforms to Structure
95.5 - 104.5 %
Clear

Colorless to Faint Yellow
≥ 95.5 %

Specification Date : 04/29/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.